

Termite Nuclei Dissociation and FACS Sorting

02/06/2023-06/06/2023

Author: Catherine Gatt, Fabilab (<https://fabilab.org>), UNSW Sydney, Australia

Experimental Overview

Goals:

1. Dissociate a decent number of nuclei from termites.
2. Remove termite debris from dissociated nuclei.
3. Achieve a concentration of nuclei appropriate for 10x Genomics sc sequencing (approx 700-1200 cells/uL)

Method:

1. Nuclei dissociation of termite nuclei using homogenisation and resuspension buffer solutions, grinding, centrifugation and straining.
2. FACS sorting using DAPI and the BD FACSAria III to collect nuclei and remove debris.

Results:

1. Current homogenisation and resuspension buffer solutions, grinding, centrifugation and straining protocol dissociates nuclei from termites.
2. FACS sorting with the BD FACSAria III collected 100K DAPI stained nuclei per run, from a sample containing 6-7 termites, and removed debris from nuclei. Resulting FACS plot looks similar to the plot of single nuclei from *Drosophila* tissues.
3. Concentration of nuclei was achieved by centrifugation.

Below I will briefly outline the protocol we used, and our results, including the FACS plots and images of nuclei under epifluorescent microscope. Please find our excel lab protocol attached to the email.

1. *Nuclei Dissociation*

Twelve *N. walkeri* soldiers and twelve workers were placed on a dissection dish and were cut into small pieces using tweezers and scissors. Approximately six termites (soldiers and workers were moved to separate tubes) where each moved into one 1.5-mL tube containing 1000 μ L of chilled homogenization buffer (Nuclease Free water, 250 mM Sucrose (nuclease free), 10 mM Tris (pH 8.0), 25 mM KCl, 5 mM MgCl₂, 0.1% Triton-x 100 and 0.1 mM DTT). The tissue was then homogenized using a pellet pestle, by twisting the pestle approximately 20 times. The nuclei were centrifuged at 1000g for 10 min at 4 °C, and the supernatant was discarded. The pellet was then resuspended in 1000 μ L of resuspension buffer (PBS and 0.5% BSA). The nuclei were then strained using 100 and 35 μ m cell strainers (see Result 1). After the previous two steps, and after FACS sorting (see Result 2), 10 μ L of the nuclei was stained with DAPI (1:1000) and visualized

under an epifluorescent microscope. Images of the cells were taken using an Iphone (1), and microscope in the DAPI (2) and brightfield (3) channels.

2. FACS Sorting

FACS sorting was performed on the BD FACSAria III. Approximately 100 μ L of nuclei suspension was run as an unstained control (see Result 3). DAPI was added (1:1000) to the remaining nuclei and was captured into a FACS tube containing 50 μ L of PBS, excluding debris and doublets, and sorting was run until (1) there was no sample left (Tube 2) or 100K events had occurred (Tube 1, 3 and 4). (see Result 3). Subsequently, the concentration of nuclei was determined with a hemocytometer and counting the number under a 20 \times objective of an epifluorescence microscope (see Result 2). The nuclei were centrifuged at 1000g for 6 min at 4 $^{\circ}$ C, and approx. 600 μ L of supernatant was discarded. The pellet was resuspended in the remaining supernatant.




Sort Settings:

Experiment : CG20230606		Sort Report		Report Date : 2023.06.06 at 16:39:24	
Specimen : Specimen_001				Device : 2 Tube	
Tube : T4-W2				User ID : catgat	
Sort Layout : Sort Layout_001				Cytometer : FACSAriaIII (1)	
Application : FACSDiva Version 8.0.1					
Sort Settings					
Sort Setup	100 micron	Precision		Purity	
Frequency	31.1	Yield Mask		32	
Amplitude	4.2	Purity Mask		32	
Phase	0.00	Phase Mask		0	
Drop Delay	28.52	Single Cell		Off	
Attenuation	Off	Plates Voltage		4,500	
Sweet Spot	On	Voltage Centering		200	
First Drop	220	Sheath Pressure		20.00	
Target Gap	10				
Side Stream Voltage (%)					
Far Left	Left	Right		Far Right	
0.00	32.00	0.00		0.00	
Neighboring Drop Charge (%)					
2nd		3rd		4th	
15.00		6.00		1.00	




Results

1. Nuclei after Straining

a. 100 μ m cell strainer





b. 35 μ m cell strainer






2. Nuclei After FACS

a. Before Concentration




b. After Concentration



3. FACS Plots

a. Unstained Control, 02/06/2023

BD FACSDiva 8.0.1



b. Tube 1 (soldiers), 06/06/2023

c. **Single nuclei from *Drosophila* tissues**

i. Source: <https://star-protocols.cell.com/protocols/1669>





Figure 2. Representative FACS plots of nuclei from adult fly heads

Left: plot of unstained nuclei. Right: plot of nuclei stained with Hoechst. The two populations of nuclei in the polygon are sorted for snRNA-seq.


d. Tube 2 (soldiers), 06/06/2023



Tube: T2-S2			
Population	#Events	%Parent	%Total
■ All Events	10,000	####	100.0
■ P1	1,131	11.3	11.3
■ P2	1,276	12.8	12.8





e. Tube 3 (workers), 06/06/2023

BD FACSDiva 8.0.1



f. Tube 4 (workers) 06/06/2023

BD FACSDiva 8.0.1



Tube: T4-W2			
Population	#Events	%Parent	%Total
■ All Events	10,000	####	100.0
■ P1	673	6.7	6.7
■ P2	709	7.1	7.1